








First clinically validated blood test to be used as an aid in the determination of menopausal status.




What is MenoCheck ?



Why MenoCheck ?



Interpretation of Results



Menopause category	Final menstrual period (FMP)	AMH concentration
Pre Menopause	>5 years from FMP	>100 pg/mL
Peri Menopause	<5 years from FMP	10-99.9 pg/mL
Post Menopause	At FMP or later	<10 pg/mL

Clinical Validation

Predict the trans-menopausal bone loss
Karlman et al, 2016

Predict the onset of vasomotor symptoms
Crawford et al, 2016

Predict the menopausal lipid profile change
El Khoudary et al, 2016

CLINICALLY VALIDATED IN LONGITUDINAL STUDIES

Predict the new onset of sleep disturbance
Joffe et al, 2016

Predict the onset of early menopause
Crawford et al, 2016

Predict the time to final menstrual period
Finkelstein et al, 2016

www.MenoCheck.com